

Fatal 4-MTA Intoxication: Development of a LC Tandem Mass Spectrometric Assay for Multiple Matrices

T. Decaestecker

E. De Letter, K. Clauwaert, M.P. Bouche, W. Lambert, J. Van Bocxlaer, M.Piette, E. Van den Eeckhout, C. Van Peteghem, and A. De Leenheer

Laboratory of Medical Biochemistry and Toxicology, Ghent University
Laboratory of Legal Medicine, Ghent University



Outline

- Introduction
- Case history
- Systematic Drug Screening
- Analysis of 4-MTA
 - Isolation of the compounds
 - Chromatography
 - Mass Spectrometry
- Validation of the method
- Distribution of 4-MTA
- Conclusions



Introduction

4-methylthioamphetamine (4-MTA)

- Street name: Flatliner
- New, non-neurotoxic serotonin-releasing agent (\leftrightarrow MDMA)
- Only **three** fatal overdose cases reported
 - ↳ toxicity range (Elliot):
 - 0.2 – 0.6 mg/L blood : moderately toxic
 - 0.6 – 1.5 mg/L blood : severely toxic
 - > 1.5 mg/L blood: fatal



Case history

- 27-year old drug dealer † after collapse and reanimation
- Examination:

External: skin lesions

Internal:

- ✓ Lungs: emphysema, **congestion** & edema
- ✓ Heart: **no** cardiac anomalies
- ✓ Brain: slight **congestion** & edema of the white matter

⇒ A wide variety of specimens was taken



Systematic Drug Screening

➤ Blood & urine

- ✓ **EMIT**: amphetamines, cannabinoids, caffeine & cotinine
- ✓ **HPLC-DAD, GC-MS, GC-NPD**: 4-MTA, MDMA & caffeine

➤ LC-MS/MS on a variety of matrices

- ↳ Detailed **distribution study** in various blood and tissue samples



Analysis of 4-MTA

Isolation of the compounds

- 1mL sample (blood, urine or tissue homogenate)
- Internal standardisation (**phentermine**)
- **L/L** extraction
 - pH 9.5 (K_2CO_3)
 - Hexane/ethyl acetate (70/30)
 - **MeOH/HCl**
 - Evaporation of organic phase





Chromatography

- Hypersil BDS phenyl column (narrow-bore)
- Gradient elution (NH_4Ac in H_2O & MeOH/AcCN)



Mass spectrometry

- Q-TOF – ESI⁺
- Parameter tuning
 - ⇒ 🖐️ **cone** (20V) & **collision energy** (11-15eV)
- Quantification: reconstructed **mass fragmentograms**





Precursor mass 182.1
Product mass 165



Precursor mass 150.1
Product masses 91 + 133.1




Method Validation Parameters (4-MTA)

	Blood	Urine	Tissue
Linearity			
Range ($\mu\text{g/L}$)	5-2500	5-2500	125-2500
Intercept	0.0049	0.0361	-0.0463
Slope	0.0030	0.0032	0.0026
R^2	0.9992	0.9992	0.9987
Total Reproducibility (n=5) CV%	<14%	<15%	<9%
LOD ($\mu\text{g/L}$)	2.5	2.5	50
LLQ ($\mu\text{g/L}$)	5	5	125
ULQ ($\mu\text{g/L}$)	2500	2500	2500



Distribution of 4-MTA (mg/L or mg/kg)



Conclusions

- Developed LC-MS/MS method
 - Selective
 - Sensitive (LOD 2.5 µg/L, except for 4-MTA in tissue homogenate)
 - Uncomplicated sample pretreatment
- Case:
 - Extensively examined
 - Interesting for cases without blood available
 - Contribute to distribution studies

